പു

Surface Wa

The surface Water Absorption Te used to determine the fabric, such as towels, bathrobes, bath mats wa beachwear, etc, to rapidly retain liquid water from s human skin, dishes and Specimens are placed at a the base of the instrume flows down on the s specimen, the amount of water retained in each specimen is measured. Six specimens are tester, three on the face side and three on the back side. The six observations are averaged to determine the surface water absorption of the fabric

ption Tester



Technical details

Specimen 180 mm \emptyset Pour spout 30 ± 2 mm

Angle of holder 60°

Waterflow 50 ml within 8 s

Test distance 2 - 10 mm

Weight 10 kg

Dimensions 36 x 26 x 65 cm

Standards

ASTM D4772